KN-TS10

## Travel scale

<b>Description</b>	<b>English</b>
The travel scale is used to measure the body weight of human beings.	

<b>Power supply</b>	
The device is powered with one lithium battery (CR2032). If the display shows "LO", replace the battery.	

- Use**
- Extract the display.
  - Place the device on a stable, flat surface.
  - To reset the device, press the platform with your foot. The display shows "0.0". The device switches on automatically. The device is ready for use.
  - Step on the platform. Place your feet evenly on the platform without moving. The display shows your weight.
  - Step off the platform. The device switches off automatically.
  - Press the unit-of-measurement button on the back of the device to adjust the unit of measurement: kilogram (kg), pound (lb) or stone (st).

- Safety**
- Read the manual carefully before use. Keep the manual for future reference.
  - Only use the device for its intended purposes. Do not use the device for other purposes than described in the manual.
  - Do not use the device if any part is damaged or defective. If the device is damaged or defective, replace the device immediately.

- Battery safety**
- Use only the batteries mentioned in the manual.
  - Do not install batteries in reverse polarity.
  - Do not short-circuit or disassemble the batteries.
  - Do not expose the batteries to water.
  - Do not expose the batteries to fire or excessive heat.
  - Batteries are prone to leakage when fully discharged. To avoid damage to the product, remove the batteries when leaving the product unattended for longer periods of time.
  - If liquid from the batteries comes into contact with skin or clothing, immediately rinse with fresh water.

## Cleaning and maintenance

- Warning!**
- Do not use cleaning solvents or abrasives.
  - Do not clean the inside of the device.
  - Do not attempt to repair the device. If the device does not operate correctly, replace it with a new device.
- Clean the outside of the device using a soft, damp cloth.**

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